## HEIZKISSEN HK 894 CB

Verwarmingskussen • Coussin chauffant • Almohadilla eléctrica • Almofada eléctrica  
Cuscino termico • Varmepute • Heating pad • Poduszka ogrzewająca  
Nahřivací dečka • Melegítőpárna • Электрическая грелка

## Bedienungsanleitung/Garantie

Gebruiksaanwijzing/Garantie • Mode d'emploi/Garantie • Instrucciones de servicio/Garantía  
Manual de instruções/Garantia • Istruzioni per l'uso/Garanzia • Bruksanvisning/Garanti  
Instruction Manual/Guarantee • Instrukcja obsługi/Gwarancja • Návod k použití/Záruka  
Használati utasítás/Garancia • Руководство по эксплуатации/Гарантия

## DEUTSCH

### Inhalt

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## Please keep for subsequent use!

### WARNING:

- This device is not suitable for use in hospitals.
- During use the user must not fall asleep.
- Do not cover up the thermostat. This may lead to malfunctions.
- If the device is used for a prolonged period on the highest setting, burning of the skin may result.
- The product must not be used for small children, frail individuals of those sensitive to heat.
- If the product is used by children, the thermostat should be set by a parent or supervisor. Use of the device must be monitored!
- Do not use the product if it is wet or in rooms with a high level of humidity, such as bathrooms etc.
- After use, do not fold the product too tightly in order not to damage internal wiring!
- The product should be inspected regularly for signs of wear or damage. If any such signs are found, or if the device has been used incorrectly, it must be returned to the dealer or service department before being used again.
- Do not use a damaged device, and disconnect it immediately from the mains if it is already in use.
- If the device does not work because a fuse has been activated or triggered, please return the product to our service department.

## General Safety Instructions

Read the operating instructions carefully before putting the appliance into operation and keep the instructions including the warranty, the receipt and, if possible, the box with the internal packing. If you give this device to other people, please also pass on the operating instructions.

- The appliance is designed exclusively for private use and for the envisaged purpose. This appliance is not fit for commercial use. Do not use it outdoors (except if it is designed to be used outdoors). Keep it away from sources of heat, direct sunlight, humidity (never dip it into any liquid) and sharp edges. Do not use the appliance with wet hands. If the appliance is humid or wet, unplug it immediately.
- When cleaning or putting it away, switch off the appliance and always pull out the plug from the socket (pull on the plug, not the cable) if the appliance is not being used and remove the attached accessories.
- Do **not** operate the machine without supervision. If you should leave the workplace, always switch the machine off or remove the plug from the socket (pull the plug itself, not the lead).

- Check the appliance and the cable for damage on a regular basis. Do not use the appliance if it is damaged.
- Do not try to repair the appliance on your own. Always contact an authorized technician. To avoid the exposure to danger, always have a faulty cable be replaced only by the manufacturer, by our customer service or by a qualified person and with a cable of the same type.
- Use only original spare parts.

## Children and Frail Individuals

- In order to ensure your children's safety, please keep all packaging (plastic bags, boxes, polystyrene etc.) out of their reach.

**Caution!** Do not allow small children to play with the foil as there is a **danger of suffocation!**

- In order to protect children or frail individuals from the hazards of electrical equipment, please ensure that this device is only used under supervision. This device is not a toy. Do **not** allow small children to play with it.

## Symbols in these Instructions for Use

Important information for your safety is specially marked. It is essential to comply with these instructions in order to avoid accidents and prevent damage to the machine:

### WARNING:

This warns you of dangers to your health and indicates possible injury risks.

### CAUTION:

This refers to possible hazards to the machine or other objects.

**NOTE:** This highlights tips and information.

## Symbols on the Product

You will also find warning and information symbols on the product:

**Read the operating instructions!**

**Do not insert any needles!**

**Do not use when folded together or not spread out!**

## Overview of the Components

- 1 Heating pad
- 2 Thermostat
- 3 Control lamp
- 4 Cover

## Instructions for Use

### Use

The product should preferably be used on the following parts of the body: shoulders, arms, legs, calves and abdomen.

#### **WARNING:**

The heating pad must always be used inside the cover (4). This ensures that a minimum distance is maintained which protects against overheating.

### Electrical Connection

- Please ensure that the device voltage (see type plate) and mains voltage are the same.
- Insert the mains plug into a properly installed 230 V, 50 Hz socket with earthing contact.

### Switching on

- Move the regulator on the thermostat (2) to the desired setting.
- The control lamp (3) lights up and shows that the device is working.

#### **NOTE:**

On the first few occasions that this device is used unusual odours may be produced. These are harmless and disappear the more the device is used.

### Thermostat Settings

You can select from the following settings:

Setting	Function
0	Off
1	Low. Slight warmth.
2	Medium
3	High

 **NOTE:** You can read off the setting on the control lamp (3).

### Switching off

- Turn the thermostat control to 0 (off).
- Only then should you remove the plug from the mains socket.

### Cleaning

#### **WARNING:**

Always remove the mains plug from the socket before cleaning the device.

### Cover

- After you have removed the heating pad you can wash the cover (4) in the washing machine.
- Please note the care instructions on the label.

### Heating Pad with Thermostat

#### **CAUTION:**

- The heating pad with thermostat must not be washed in the washing machine or dry-cleaned.
- Do not immerse the components in water.
- The thermostat must not be allowed to become wet.
- Do not use any abrasive cleaners or solvents.
- Wipe the heating pad off with a damp cloth. In the case of stubborn stains you can use a gentle solution of soapy water.
- Dry off the heating pad after cleaning.
- You should only wipe off thermostat with a dry cloth.

### Storage

- Remove the mains plug from the socket and allow the heating pad to cool down for approx. 20 minutes.
- Remove the cover. You can now wash it if necessary. Please note the care instructions on the label.
- Fold the heating pad together gently and store the entire product in the storage bag.

#### **CAUTION:**

- Store the product in a dry place.
- Do not place any heavy objects on the device.

## Technical Data

Model:..... HK 894 CB

Power supply:..... 230 V, 50 Hz

Power consumption:..... 100 W

Protection class:..... II

This device has been tested according to all relevant current CE guidelines, such as electromagnetic compatibility and low voltage directives, and has been constructed in accordance with the latest safety regulations.

Subject to technical changes without prior notice!

## Guarantee

The device supplied by our Company is covered by a 24 month guarantee starting on the date of purchase (receipt).

During the life of the guarantee any fault of the device or its accessories ascribable to material or manufacturing defects will be eliminated free of charge by repairing or, at our discretion, by replacing it. The guarantee services do not entail an extension of the life of the guarantee nor do they give rise to any right to a new guarantee!

Proof of the guarantee is provided by the proof of purchase. Without proof of purchase no free replacement or repair will be carried out.

If you wish to make a claim under the guarantee please return the **entire** machine in the **original packaging** to your dealer together with the receipt.

Defects to consumables or parts subject to wearing as well as cleaning, maintenance or the replacement of said parts are not covered by the guarantee and hence are to be paid!

The guarantee lapses in case of unauthorized tampering.

### **After the expiry of the guarantee**

After the expiry of the guarantee repairs can be carried by the competent dealer or repair service against the payment of the ensuing costs.

### **Meaning of the “Dustbin” Symbol**

Protect our environment: do not dispose of electrical equipment in the domestic waste.

Please return any electrical equipment that you will no longer use to the collection points provided for their disposal.

This helps avoid the potential effects of incorrect disposal on the environment and human health.

This will contribute to the recycling and other forms of reutilisation of electrical and electronic equipment.

Information concerning where the equipment can be disposed of can be obtained from your local authority.

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C. Bomann GmbH